

Needle-Knife Sphincterotome disposable




STERILE EO

disposable Needle-Knife Sphincterotome

- flexible rounded needle 0.9 mm Ø
- exposure 6 mm
- tapered tip
- colored marked
- lateral Luer Lock connector for contrast injektion

flexible Needle



Order No.	Tube diameter	Needle diameter	Needle Stop	Effektive length	Packing unit
OE11018N6-09	1.8 mm	0.9 mm	6 mm	200 cm	5 pieces

Subject to change.

latest version: 10/2015