

pre-loaded Stent-Set



STERILE EO



pre-loaded biliary PE-Stent 10 Fr. ;
available in 5 Stent Lengths

disposable Stent-Set

contains:

- Stent Application Set
- pre-loaded PE-Stent (straight)
- Luer-Lock-Connector

Order No.	Pusher Ø	Guiding catheter Ø	Total length	PE-Stent	Packing unit
SET-01-10-05	10 Fr.	6 Fr.	220 cm	10 Fr. / 5 cm	1 piece
SET-01-10-07	10 Fr.	6 Fr.	220 cm	10 Fr. / 7 cm	1 piece
SET-01-10-09	10 Fr.	6 Fr.	220 cm	10 Fr. / 9 cm	1 piece
SET-01-10-11	10 Fr.	6 Fr.	220 cm	10 Fr. / 11 cm	1 piece
SET-01-10-13	10 Fr.	6 Fr.	220 cm	10 Fr. / 13 cm	1 piece

Änderungen vorbehalten.

latest version: 01/2016