


Single-use	Cytology Brushes	5.1
Single-use	Cleaning Brushes	5.2 - 5.3
Reusable	Cleaning Brushes	5.4



Ø 2.6
mm



STERILE EO

Cytology Brush

- double lumen biliary cytology brush
- for guide wire .035 inch



single-use

Cytology Brushes

length

packing unit

Order No.

180 cm

5 pieces

E 42 22 DL-C



Ø 1.8
mm



STERILE EO

Cytology Brush

- metal ball for closing the tube



length

packing unit

Order No.

120 cm


10 pieces

NE 42 18-A

180 cm

10 pieces

NE 42 18-C



Ø 2.3
mm



STERILE EO

Cytology Brush

- metal ball for closing the tube



length

packing unit

Order No.

230 cm

10 pieces

NE 42 22-G



Port Cleaning Brush

- double sided
- for cleaning of endoscope valves



small brush Ø	large brush Ø	packing unit	Order No.
5.0 mm	12.0 mm	25 pieces	ZB0050


Ø 2.2 mm



Cleaning Brush Set

- Brush 1: working channels
- Brush 2: valve / Port Brush

set components	total length	brush diameter	tube diameter	packing unit	Order No.
Brush 1	230 cm	7 mm	1.8 mm	25 pieces	NEBS 41221-G-SET01
Brush 2		5 and 12 mm	-		
Brush 1	180 cm	5 mm	1.8 mm	25 pieces	NEBS 41221-C-SET01
Brush 2		5 and 12 mm	-		



Ø 1.8
mm



Cleaning Brush

- double sided
- for cleaning of endoscope working channels



length	packing unit	length of brushes	brush Ø 5.0 mm
120 cm	10 pieces	20 mm	EBS 41 181-A

Ø 2.2
mm



Cleaning Brush

- double sided
- for cleaning of endoscope working channels



length	packing unit	length of brushes	brush Ø 5.0 mm
180 cm	10 pieces	20 mm	EBS 41 221-C
230 cm	10 pieces	20 mm	EBS 41 221-G

Cleaning Brush

- double sided
- for cleaning of endoscope working channels
- perfectly rounded brush heads



length	packing unit	length of brushes	Colour	brush Ø 5.0 mm
180 cm	25 pieces	20 mm	yellow	NEBS 41 221-C
230 cm	25 pieces	20 mm	green	NEBS 41 221-G
280 cm	25 pieces	20 mm	lilac	NEBS 41 221-M

Ø 1.8
mm

134 °C
273 °F

single / double sided

- for cleaning of endoscope working channels



brush Ø 3.0 mm

length	single sided, incl. handle	double sided
120 cm	41 181-A	BS 41 181-A
160 cm	41 181-B	BS 41 181-B
180 cm	41 181-C	BS 41 181-C
280 cm	41 181-M	

Ø 2.2
mm

134 °C
273 °F

single / double sided

- for cleaning of endoscope working channels



brush Ø 5.0 mm

length	single sided, incl. handle	double sided
120 cm	41 221-A	BS 41 221-A
160 cm	41 221-B	BS 41 221-B
180 cm	41 221-C	BS 41 221-C
230 cm	41 221-G	BS 41 221-G
280 cm	41 221-M	

single / double sided

- for cleaning of endoscope working channels



brush Ø 7.0 mm

length	single sided, incl. handle	double sided
120 cm	41 401-A	BS 41 401-A
160 cm	41 401-B	BS 41 401-B
180 cm	41 401-C	BS 41 401-C
230 cm	41 401-G	BS 41 401-G