

Single-use	Injection Needles, PRO-SERIES	4.1 - 4.2
Single-use	Injection Needles, STANDARD-SERIES	4.3
Single-use	Transbronchial Aspiration Needles	4.4
Single-use	Fibrin Application Needles	4.5

Injection Needles





- with metal tip for safe needle exposure

Ø 2.3 mm



STERILE EO



Ø 0.7 mm

length	packing unit	4 mm	5 mm	6 mm
160 cm	5 pieces	NET 25 22-B4		NET 25 22-B6
180 cm	5 pieces		NET 25 22-C5	
230 cm	5 pieces	NET 25 22-G4	NET 25 22-G5	NET 25 22-G6



Ø 1.0 mm

length	packing unit	4 mm	6 mm
160 cm	5 pieces	NET 26 22-B4	NET 26 22-B6
230 cm	5 pieces	NET 26 22-G4	NET 26 22-G6

single-use

Injection Needles



- with metal tip for safe needle exposure

Ø 2.3 mm



STERILE EO



Ø 0.5 mm

length	packing unit	4 mm	6 mm
230 cm	5 pieces	NET 24 22-G4	NET 24 22-G6

Ø 1.8 mm



STERILE EO



Ø 0.5 mm

length	packing unit	4 mm	6 mm
180 cm	5 pieces	NET 25 18-C4	NET 25 18-C6

Ø 2.3 mm



STERILE EO



Ø 0.7 mm

length	packing unit	4 mm	6 mm
160 cm	5 pieces		ET 25 22-B6
180 cm	5 pieces	ET 25 22-C4	ET 25 22-C6
230 cm	5 pieces	ET 25 22-G4	ET 25 22-G6
280 cm	5 pieces		ET 25 22-M6
300 cm	5 pieces		ET 25 22-O6



Ø 0.9 mm

length	packing unit	6 mm
160 cm	5 pieces	ET 26 22-B6
180 cm	5 pieces	ET 26 22-C6
230 cm	5 pieces	ET 26 22-G6

Ø 1.8 mm



STERILE EO



Ø 0.6 mm

length	packing unit	4 mm	6 mm
180 cm	5 pieces	ET 25 18-C4	ET 25 18-C6
280 cm	5 pieces		ET 25 18-M6


Ø 1.9 mm



STERILE EO

Transbronchial Aspiration Needle

- Outer **metal sheath** for enhanced punction force



Ø 0.8 mm

Ø 0.8 mm

length	packing unit	12 mm	15 mm
120 cm	5 pieces	MEB 26 18-A12	MEB 26 18-A15


Ø 2.3 mm



STERILE EO

Transbronchial Aspiration Needle

- Needle with lateral opening, Handle incl. needle safety device



Ø 0.7 mm

Ø 0.9 mm

length	packing unit	15 mm	15 mm
120 cm	5 pieces	EB 25 22-A15	EB 26 22-A15


Ø 1.8 mm



STERILE EO

Transbronchial Aspiration Needle

- Needle with lateral opening, Handle incl. needle safety device



Ø 0.7 mm

Ø 0.9 mm

length	packing unit	15 mm	12 mm
120 cm	5 pieces	EB 25 19-A15	EB 26 18-A12

Ø 2.6 mm



STERILE EO

Fibrin Application Needle

- double-lumen
- single needle

length	6 mm
160 cm	ET 27 2620-B
180 cm	ET 27 2620-C
230 cm	ET 27 2620-G

Ø 0.9 mm



Fibrin Application Needle

- double-lumen
- coaxial needle

length	6 mm
160 cm	ET 27 2630-B
180 cm	ET 27 2630-C
230 cm	ET 27 2630-G

Ø 1.2 mm



single-use

Fibrin Application Needles