

Single-use	Polypectomy Snares, rotatable	3.1 - 3.2
Reusable	Polypectomy Snares, fixed handle	3.3
Reusable	Polypectomy Snares, rotatable	3.4
Reusable	Special shapes	3.5

Polypectomy Snares



Ø 2.3
mm



STERILE EO



rotatable

oval loop

- braided wire



length	packing unit	15 mm	25 mm	35 mm
180 cm	5 pieces	NOE 34 2214-C	NOE 34 2216-C	NOE 34 2217-C
230 cm	5 pieces	NOE 34 2214-G	NOE 34 2216-G	NOE 34 2217-G

asymmetrical loop

- braided wire



length	packing unit	15 mm	25 mm
230 cm	5 pieces	NOE 3222-G	NOE 3322-G

hexagonal loop

- braided wire



length	packing unit	15 mm	25 mm	35 mm
230 cm	5 pieces	NOE 372214-G	NOE 372216-G	NOE 372217-G

oval loop

- monofilament wire



length	packing unit	6 mm	10 mm
230 cm	5 pieces	NOE 34 2212M-G	NOE 34 2213M-G

length	packing unit	15 mm	25 mm	35 mm
180 cm	5 pieces	NOE 34 2214M-C	NOE 34 2216M-C	NOE 34 2217M-C
230 cm	5 pieces	NOE 34 2214M-G	NOE 34 2216M-G	NOE 34 2217M-G



single-use

Polypectomy Snares

Ø 1.8
mm



STERILE EO



rotatable

oval loop

- Pediatric Snare
- braided wire



length	packing unit	25 mm
180 cm	5 pieces	NOE 34 1816-C



recommended equipment

Active Cords

HF generator	Order No.
ERBE E- and T-Series	640 300
ERBE ICC	640 500
OLYMPUS HF 120/300	640 500



Ø 2.3 mm

134 °C
273 °F

fixed handle

asymmetrical loop

- braided wire



	length	15 mm	25 mm	15 mm	25 mm
180 cm		OK 30 22-C	OK 31 22-C	OK 32 22-C	OK 33 22-C
230 cm		OK 30 22-G	OK 31 22-G	OK 32 22-G	OK 33 22-G

oval loop

- braided wire



	length	20 mm	25 mm	35 mm	20 mm	25 mm	35 mm
180 cm		OK 34 2205-C	OK 34 2206-C	OK 34 2207-C	OK 34 2215-C	OK 34 2216-C	OK 34 2217-C
230 cm		OK 34 2205-G	OK 34 2206-G	OK 34 2207-G	OK 34 2215-G	OK 34 2216-G	OK 34 2217-G

hexagonal

- braided wire



	length	20 mm	25 mm	35 mm	20 mm	25 mm	35 mm
180 cm		OK 38 2205-C	OK 38 2206-C	OK 38 2207-C	OK 37 2205-C	OK 37 2206-C	OK 37 2207-C
230 cm		OK 38 2205-G	OK 38 2206-G	OK 38 2207-G	OK 37 2205-G	OK 37 2206-G	OK 37 2207-G

oval loop

- monofilament wire



	length	20 mm	25 mm	35 mm
180 cm		OK 34 2215M-C	OK 34 2216M-C	OK 34 2217M-C
230 cm		OK 34 2215M-G	OK 34 2216M-G	OK 34 2217M-G



Ø 2.3 mm

134 °C
273 °F



without handle, rotatable

asymmetrical loop

- braided wire

insulated tip



standard tip



length	15 mm	25 mm	15 mm	25 mm
180 cm	D 30 22-C	D 31 22-C	D 32 22-C	D 33 22-C
230 cm	D 30 22-G	D 31 22-G	D 32 22-G	D 33 22-G

oval loop

- braided wire

insulated tip



standard tip



length	20 mm	25 mm	35 mm	20 mm	25 mm	35 mm
180 cm	D 34 2205-C	D 34 2206-C	D 34 2207-C	D 34 2215-C	D 34 2216-C	D 34 2217-C
230 cm	D 34 2205-G	D 34 2206-G	D 34 2207-G	D 34 2215-G	D 34 2216-G	D 34 2217-G

hexagonal

- braided wire

insulated tip



standard tip



length	20 mm	25 mm	35 mm	20 mm	25 mm	35 mm
180 cm	D 38 2205-C	D 38 2206-C	D 38 2207-C	D 37 2205-C	D 37 2206-C	D 37 2207-C
230 cm	D 38 2205-G	D 38 2206-G	D 38 2207-G	D 37 2205-G	D 37 2206-G	D 37 2207-G

oval loop

- monofilament wire

standard tip



length	20 mm	25 mm	35 mm
180 cm	D 34 2215M-C	D 34 2216M-C	D 34 2217M-C
230 cm	D 34 2215M-G	D 34 2216M-G	D 34 2217M-G

to be used with:



Order No. 620 300

Ø 2.3 mm

134 °C
273 °F

incl. handle slider

asymmetrical loop, extra large

- braided wire



insulated tip

standard tip

length	40 mm	60 mm	40 mm	60 mm
180 cm	OK 65 22401-C	OK 65 22601-C	OK 65 2240-C	OK 65 2260-C
230 cm	OK 65 22401-G	OK 65 22601-G	OK 65 2240-G	OK 65 2260-G



fixed handle

special type

- thermo-coated

length	35 mm
230 cm	OK 34 2237-G



fixed handle

special type

- semi-rigid

length	20 – 35 mm
230 cm	OK 34 2267-G



recommended equipment

Active Cords

HF generator	Order No.
ERBE E- and T-Series	640 300
ERBE ICC	640 500
OLYMPUS HF 120/300	640 500



universal handle

	Order No.
(for snares, which are supplied without handle)	620 300

